




BIOMEDORTHO
define excellence



LOCKING

PLATE SYSTEM

Ø 2.4 MM / 2.7 MM VA - TWO COLUMN VOLAR DISTAL RADIUS PLATE



Ø 2.4 MM LCP DISTAL RADIUS PLATE STRAIGHT - DORSAL



Ø 2.4 MM LCP T DISTAL RADIUS PLATE RIGHT ANGLE - DORSAL



Ø 2.4 MM LCP L DISTAL RADIUS PLATE - DORSAL



Ø 2.4 MM LCP L DISTAL RADIUS PLATE - DORSAL



Ø 2.4 MM LCP L DISTAL RADIUS PLATE OBLIQUE ANGLE - DORSAL



**Ø 2.4 MM LCP DISTAL RADIUS
VOLAR PLATE - VOLAR**



**Ø 2.4 MM LCP DISTAL RADIUS
EXTRA ARTICULAR PLATE - VOLAR**



**Ø 2.4 MM LCP DISTAL RADIUS
EXTRA ARTICULAR PLATE - VOLAR**



**Ø 2.4 MM LCP DISTAL RADIUS
BUTTRESS PLATE - VOLAR**



**Ø 2.4 MM LCP DISTAL RADIUS
COLUMN VOLAR PLATE - VOLAR**




**Ø 2.4 MM LCP DISTAL RADIUS
COLUMN VOLAR PLATE - VOLAR**




Ø 2.4 MM LCP RADIAL HEAD RIM PLATE



Ø 2.4 MM LCP RADIAL HEAD NECK PLATE



Ø 2.4 MM LOCKING SCREW - SELF TAPPING



Ø 2.7 MM LOCKING SCREW - SELF TAPPING




Ø 2.7 MM CORTICAL SCREW - SELF TAPPING



3.5 MM

Ø 3.5 MM LCP SUPERIOR ANTERIOR CLAVICLE PLATE WITHOUT LATERAL EXTENSION



Ø 3.5 MM LCP CLAVICLE HOOK PLATE




Ø 3.5 MM LCP PROXIMAL HUMERUS PLATE (LONG PHILOS)



Ø 3.5 MM LCP SUPERIOR ANTERIOR CLAVICLE PLATE WITH LATERAL EXTENSION



Ø 3.5 MM LCP PROXIMAL HUMERUS PLATE (SHORT PHILOS)



Ø 3.5 MM LCP PERIARTICULAR PROXIMAL HUMERUS PLATE



Ø 2.7 / 3.5 MM LCP POSTERO -LATERAL DISTAL HUMERUS PLATE WITH SUPPORT

SUPPORTED SCREW

3.5 MM LOCKING SCREW - SELF TAPPING
2.7 MM LOCKING SCREW - SELF TAPPING
3.5 MM LCP CANCELLOUS SCREW - FULL THREADED
4.0 MM CANCELLOUS SCREW - FULL THREADED
3.5 MM CORTICAL SCREW - SELF TAPPING

➤ **HOLES** 3 to 14
➤ **TYPE** SS & TT
➤ **ANGLED** LEFT & RIGHT
➤ **WIDTH** 11.0 MM
➤ **THICKNESS** 2.5 MM
➤ **LENGTH** 116 MM
FOR 7 HOLES

3.5 MM LOCKING SCREW

2.7 MM LOCKING SCREW

3.5 MM LCP CANCELLOUS SCREW

3.5 MM CORTICAL SCREW

4.0 MM CANCELLOUS SCREW

Ø 3.5 MM LCP MEDIAL DISTAL HUMERUS PLATE

SUPPORTED SCREW

3.5 MM LOCKING SCREW - SELF TAPPING
2.7 MM LOCKING SCREW - SELF TAPPING
3.5 MM LCP CANCELLOUS SCREW - FULL THREADED
4.0 MM CANCELLOUS SCREW - FULL THREADED
3.5 MM CORTICAL SCREW - SELF TAPPING

➤ **HOLES** 3 to 14
➤ **TYPE** SS & TT
➤ **ANGLED** LEFT & RIGHT
➤ **WIDTH** 11.0 MM
➤ **THICKNESS** 2.5 MM
➤ **LENGTH** 110 MM
FOR 7 HOLES

3.5 MM LOCKING SCREW

2.7 MM LOCKING SCREW

3.5 MM LCP CANCELLOUS SCREW

3.5 MM CORTICAL SCREW

4.0 MM CANCELLOUS SCREW

Ø 3.5 MM LCP SMALL PLATE

SUPPORTED SCREW

3.5 MM LOCKING SCREW - SELF TAPPING
3.5 MM LCP CANCELLOUS SCREW - FULL THREADED
4.0 MM CANCELLOUS SCREW - FULL THREADED
3.5 MM CORTICAL SCREW - SELF TAPPING

➤ **HOLES** 4 to 14
➤ **TYPE** SS & TT
➤ **WIDTH** 11.0 MM
➤ **THICKNESS** 3.4 MM
➤ **LENGTH** 137 MM
FOR 10 HOLES

3.5 MM LOCKING SCREW

3.5 MM LCP CANCELLOUS SCREW

4.0 MM CANCELLOUS SCREW

3.5 MM CORTICAL SCREW

Ø 2.7 / 3.5 MM LCP POSTERO -LATERAL DISTAL HUMERUS PLATE WITHOUT SUPPORT

SUPPORTED SCREW

3.5 MM LOCKING SCREW - SELF TAPPING
2.7 MM LOCKING SCREW - SELF TAPPING
3.5 MM LCP CANCELLOUS SCREW - FULL THREADED
4.0 MM CANCELLOUS SCREW - FULL THREADED
3.5 MM CORTICAL SCREW - SELF TAPPING

➤ **HOLES** 3 to 14
➤ **TYPE** SS & TT
➤ **ANGLED** LEFT & RIGHT
➤ **WIDTH** 11.0 MM
➤ **THICKNESS** 2.5 MM
➤ **LENGTH** 116 MM
FOR 7 HOLES

3.5 MM LOCKING SCREW

2.7 MM LOCKING SCREW

3.5 MM LCP CANCELLOUS SCREW

3.5 MM CORTICAL SCREW

4.0 MM CANCELLOUS SCREW

Ø 3.5 MM LCP EXTRA ARTICULAR DISTAL HUMERUS PLATE

SUPPORTED SCREW

3.5 MM LOCKING SCREW - SELF TAPPING
3.5 MM LCP CANCELLOUS SCREW - FULL THREADED
4.0 MM CANCELLOUS SCREW - FULL THREADED
3.5 MM CORTICAL SCREW - SELF TAPPING

➤ **HOLES** 4 to 14
➤ **TYPE** SS & TT
➤ **ANGLED** LEFT & RIGHT
➤ **WIDTH** 12.5 MM
➤ **THICKNESS** 3.8 MM
➤ **LENGTH** 194 MM
FOR 8 HOLES

3.5 MM LOCKING SCREW

3.5 MM LCP CANCELLOUS SCREW

4.0 MM CANCELLOUS SCREW

3.5 MM CORTICAL SCREW

Ø 3.5 MM LCP METAPHYSEAL DISTAL MEDIAL HUMERUS PLATE

SUPPORTED SCREW

3.5 MM LOCKING SCREW - SELF TAPPING
3.5 MM LCP CANCELLOUS SCREW - FULL THREADED
4.0 MM CANCELLOUS SCREW - FULL THREADED
3.5 MM CORTICAL SCREW - SELF TAPPING

➤ **HOLES** 3 to 13
➤ **TYPE** SS & TT
➤ **WIDTH** 11.0 MM
➤ **THICKNESS** 2.5 MM
➤ **LENGTH** 158 MM
FOR 13 HOLES

3.5 MM LOCKING SCREW

3.5 MM LCP CANCELLOUS SCREW

4.0 MM CANCELLOUS SCREW

3.5 MM CORTICAL SCREW

3.5 MM

Ø 3.5 MM LCP CLOVERLEAF PLATE



Ø 3.5 MM LCP RECONSTRUCTION PLATE (STRAIGHT)



Ø 3.5 MM LCP 1/3 TUBULAR PLATE



Ø 3.5 MM LCP OLECRANON PLATE



Ø 3.5 MM LCP T-PLATE RIGHT ANGLE



Ø 3.5 MM LCP T-PLATE OBLIQUE ANGLE



Ø 3.5 MM LCP PROXIMAL TIBIA PLATE - RAFT



Ø 3.5 MM LCP MEDIAL PROXIMAL TIBIA PLATE



Ø 3.5 MM LCP POSTERIOR MEDIAL PROXIMAL TIBIA PLATE



Ø 3.5 MM LCP LOW BEND MEDIAL DISTAL TIBIA PLATE



Ø 3.5 MM LCP MEDIAL DISTAL TIBIA PLATE WITH TAB



Ø 3.5 MM LCP ANTEROLATERAL DISTAL TIBIA PLATE




Ø 2.7 / 3.5 MM LCP LATERAL DISTAL FIBULA PLATE



Ø 3.5 MM LCP SMALL HOOK PLATE (FEBULA, TIBIA, RADIUS & ULNA)



Ø 3.5 MM LCP CALCANEAL PLATE



Ø 2.7 / 3.5 MM LCP PERI - LOCK VOLAR DISTAL RADIUS PLATE



Ø 3.5 MM LCP RECONSTRUCTION PLATE - ROUND HOLES



Ø 3.5 MM LOCKING SCREW

- SEFT TAPPING



- LENGTH 10 MM to 70 MM
- TYPE SS & TT
- DRILL BIT 2.8 MM
- DRILL SLEEVE 3.5 MM
- HEAD DIA 5.0 MM
- HEXAGONAL SOCKET 2.5 MM

Ø 2.7 MM LOCKING SCREW

- SEFT TAPPING



- LENGTH 8 MM to 40 MM
- TYPE SS & TT
- DRILL BIT 2.0 MM
- DRILL SLEEVE 2.7 MM
- HEAD DIA 3.5 MM
- HEXAGONAL SOCKET 2.0 MM

Ø 3.5 MM LOCKING CANCELLOUS SCREW - FULLY THREADED



- LENGTH 10 MM to 70 MM
- TYPE SS & TT
- STEP DRILL BIT 2.0 MM
- DRILL SLEEVE 2.7 MM
- HEAD DIA 5.0 MM
- HEXAGONAL SOCKET 2.5 MM

Ø 4.0 MM CANCELLOUS SCREW - FULLY THREADED



- LENGTH 10 MM to 70 MM
- TYPE SS & TT
- DRILL BIT & TAP SLEEVE 2.0 MM & 4.0 MM
- BONE TAP 4.0 MM
- HEAD DIA 6.0 MM
- HEXAGONAL SOCKET 2.5 MM

Ø 3.5 MM CORTICAL SCREW - SELF TAPPING

- SELF TAPPING



- LENGTH 10 MM to 70 MM
- TYPE SS & TT
- DRILL BIT 2.5 MM
- DRILL BIT & TAP SLEEVE 2.5 MM & 3.5 MM
- BONE TAP 3.5 MM
- HEAD DIA 6.0 MM
- HEXAGONAL SOCKET 2.5 MM

Ø 2.7 MM CORTICAL SCREW - SELF TAPPING

- SELF TAPPING



- LENGTH 8 MM to 40 MM
- TYPE SS & TT
- DRILL BIT 2.0 MM
- DRILL BIT & TAP SLEEVE 2.0 MM & 2.7 MM
- BONE TAP 2.7 MM
- HEAD DIA 4.0 MM
- HEXAGONAL SOCKET 2.0 MM

5.0 MM

Ø 4.5 / 5.0 / 6.5 MM LCP CONDYLAR BUTTRESS PLATE



Ø 4.5 / 5.0 / 6.5 MM LCP PROXIMAL FEMUR HOOK PLATE



Ø 4.5 / 5.0 MM LCP LATERAL DISTAL FEMUR PLATE - TOMOFIX



Ø 4.5 / 5.0 / 6.5 MM LCP PROXIMAL FEMUR PLATE (TROCHANTERIC PLATE)



Ø 4.5 / 5.0 MM LCP DISTAL FEMUR PLATE




Ø 4.5 / 5.0 MM LCP BROAD PLATE



Ø 4.5 / 5.0 MM LCP PROXIMAL LATERAL TIBIA PLATE - HOCKEY PLATE



Ø 4.5 / 5.0 MM LCP LATERAL HIGH TIBIA PLATE - TOMOFIX



Ø 4.5 / 5.0 MM LCP MEDIAL HIGH TIBIAL OSTEOTOMY PLATE - TOMOFIX



Ø 4.5 / 5.0 MM LCP T PLATE



Ø 4.5 / 5.0 MM LCP RECONSTRUCTION PLATE (STRAIGHT)



Ø 3.5 / 4.5 / 5.0 MM LCP METAPHYSEAL PLATE



5.0 MM

Ø 4.5 / 5.0 MM LCP T-BUTTRESS PLATE



Ø 4.5 / 5.0 MM LCP L-BUTTRESS PLATE



Ø 4.5 / 5.0 MM LCP NARROW PLATE



Ø 4.5 / 5.0 MM 135° LCP DHS PLATE LONG BARREL



Ø 4.5 / 5.0 MM 130° LCP DHS PLATE LONG BARREL



Ø 4.5 / 5.0 MM 135° LCP DHS PLATE SHORT BARREL





Ø 4.5 / 5.0 MM 130° LCP DHS PLATE SHORT BARREL



Ø 5.0 MM LOCKING SCREW - SELF TAPPING



Ø 3.5 MM LOCKING SCREW - SELF TAPPING



Ø 4.5 / 5.0 / 6.5 MM 95° LCP DCS PLATE



Ø 5.0 MM LOCKING CANCELLOUS SCREW - FULLY THREADED



Ø 3.5 MM LOCKING CANCELLOUS SCREW - FULLY THREADED



Ø 4.0 MM CANCELLOUS SCREW


- FULLY THREADED



- LENGTH 10 MM to 70 MM
- TYPE SS & TT
- DRILL BIT & TAP SLEEVE 2.0 MM & 4.0 MM
- BONE TAP 4.0 MM
- HEAD DIA 6.0 MM
- HEXAGONAL SOCKET 2.5 MM


Ø 6.5 MM CANCELLOUS SCREW

- FULLY THREADED



- LENGTH 25 MM to 110 MM
- TYPE SS & TT
- DRILL BIT 3.2 & 4.5 MM
- DRILL BIT & TAP SLEEVE 3.2 MM & 6.5 MM
- BONE TAP 6.5 MM
- HEAD DIA 8.0 MM
- HEXAGONAL SOCKET 3.5 MM


Ø 6.5 MM LOCKING CANCELLOUS SCREW - FULLY THREADED



- LENGTH 25 MM to 110 MM
- TYPE SS & TT
- STEP DRILL BIT 3.5 & 4.3 MM
- DRILL SLEEVE 6.5 MM
- HEAD DIA 8.0 MM
- HEXAGONAL SOCKET 3.5 MM

Ø 3.5 MM CORTICAL SCREW


- SELF TAPPING



- LENGTH 10 MM to 70 MM
- TYPE SS & TT
- DRILL BIT 2.5 MM
- DRILL BIT & TAP SLEEVE 2.5 MM & 3.5 MM
- BONE TAP 3.5 MM
- HEAD DIA 6.0 MM
- HEXAGONAL SOCKET 2.5 MM

Ø 4.5 MM CORTICAL SCREW

- SELF TAPPING




- LENGTH 14 MM to 70 MM
- TYPE SS & TT
- DRILL BIT 3.2 MM
- DRILL BIT & TAP SLEEVE 3.2 MM & 4.5 MM
- BONE TAP 4.5 MM
- HEAD DIA 8.0 MM
- HEXAGONAL SOCKET 3.5 MM

Ø 12.5 MM DHS / DCS COMPRESSION LAG SCREW



- LENGTH 50 MM to 120 MM
- TYPE SS & TT

2.4 / 2.7 MM LOCKING - INSTRUMENT SET





2.4 / 2.7 MM LOCKING - INSTRUMENT SET

REF. NO	COMMODITY NAME	QTY
➤ BO1000INS2401	Ø 1.8 MM DRILL BIT WITH QUICK COUPLING X 5" FOR 2.4 MM LOCKING SCREW	2
➤ BO1000INS2402	Ø 2.0 MM DRILL BIT WITH QUICK COUPLING X 5" FOR 2.7 MM LOCKING & CORTICAL SCREW	2
➤ BO1000INS2403	Ø 1.8 MM / 2.4 MM LOCKING DRILL SLEEVE FOR 2.4 MM LOCKING SCREW	2
➤ BO1000INS2404	Ø 2.0 / 2.7 MM LOCKING DRILL SLEEVE FOR 2.7 MM LOCKING SCREW	2
➤ BO1000INS2405	Ø 2.0 / 2.7 MM DRILL & TAP SLEEVE COMBINED FOR 2.7 MM CORTICAL SCREW	1
➤ BO1000INS2406	Ø 2.7 MM BONE TAP WITH QUICK COUPLING FOR 2.7 MM CORTICAL SCREW	1
➤ BO1000INS2407	T-HANDLE WITH QUICK COUPLING	1
➤ BO1000INS2408	Ø 2.4 / 2.7 MM SCREW DRIVER WITH FIBRE HANDLE	1
➤ BO1000INS2409	Ø 1.5 MM THREADED GUIDE WIRE 150 MM	2
➤ BO1000INS2410	Ø 2.4 / 2.7 MM DEPTH GAUGE	1
➤ BO1000INS2411	Ø 2.4 / 2.7 MM TORQUE LIMITING SCREW DRIVER	1
➤ BO1000INS2412	Ø 2.4 / 2.7 MM LOCKING INSTRUMENTS BOX WITH IMPLANT & SCREW TRAY	1

2.7 / 3.5 MM LOCKING - INSTRUMENT SET





2.7 / 3.5 MM LOCKING - INSTRUMENT SET

REF. NO	COMMODITY NAME	QTY
➤ BO1001INS3501	Ø 2.0 MM DRILL BIT WITH QUICK COUPLING X 5" FOR 2.7 MM LOCKING & 4.0 MM CANCELLOUS SCREW	2
➤ BO1001INS3502	Ø 2.5 MM DRILL BIT WITH QUICK COUPLING X 5" FOR 3.5 MM CORTICAL SCREW	2
➤ BO1001INS3503	Ø 2.5 / 2.8 MM STEP DRILL BIT WITH QUICK COUPLING X 8" FOR 3.5 MM LOCKING CANCELLOUS SCREW	2
➤ BO1001INS3504	Ø 2.8 MM DRILL BIT WITH QUICK COUPLING X 8" FOR 3.5 MM LOCKING SCREW	2
➤ BO1001INS3505	Ø 2.0 / 2.7 MM LOCKING DRILL SLEEVE FOR 2.7 MM LOCKING SCREW	2
➤ BO1001INS3506	Ø 2.5 / 2.8 / 3.5 MM LOCKING DRILL SLEEVE FOR 3.5 LOCKING CANCELLOUS SCREW	2
➤ BO1001INS3507	Ø 2.5 / 3.5 MM DRILL & TAP SLEEVE COMBINED FOR 3.5 MM CORTICAL SCREW	1
➤ BO1001INS3508	Ø 2.0 / 4.0 MM DRILL & TAP SLEEVE COMBINED FOR 4.0 MM CANCELLOUS SCREW	1
➤ BO1001INS3509	Ø 3.5 MM BONE TAP WITH QUICK COUPLING FOR 3.5 MM CORTICAL SCREW	1
➤ BO1001INS3510	Ø 4.0 MM BONE TAP WITH QUICK COUPLING FOR 4.0 MM CANCELLOUS SCREW	1

2.7 / 3.5 MM LOCKING - INSTRUMENT SET



REF. NO	COMMODITY NAME	QTY
➤ BO1001INS3511	T-HANDLE WITH QUICK COUPLING	1
➤ BO1001INS3512	Ø 2.7 MM SCREW DRIVER WITH FIBRE HANDLE	1
➤ BO1001INS3513	Ø 3.5 MM SCREW DRIVER WITH FIBRE HANDLE	1
➤ BO1001INS3514	Ø 1.5 MM THREADED GUIDE WIRE X 150 MM	2
➤ BO1001INS3515	Ø 2.0 MM THREADED GUIDE WIRE X 150 MM	2
➤ BO1001INS3516	Ø 2.7 MM DEPTH GAUGE	1
➤ BO1001INS3517	Ø 3.5 MM DEPTH GAUGE	1
➤ BO1001INS3518	Ø 2.7 MM TORQUE LIMITING SCREW DRIVER	1
➤ BO1001INS3519	Ø 3.5 MM TORQUE LIMITING SCREW DRIVER	1
➤ BO1001INS3520	Ø 2.7 / 3.5 MM LOCKING INSTRUMENTS BOX WITH IMPLANT & SCREW TRAY	1



BIOMEDORTHO
define excellence

**4.5 / 5.0 mm LCP
Instrument set**

4.5 / 5.0 MM LOCKING - INSTRUMENT SET



P-20



www.biomedortho.com

BIOMEDORTHO

4.5 / 5.0 MM LOCKING - INSTRUMENT SET

REF. NO	COMMODITY NAME	QTY
➤ BO1002INS5001	Ø 3.2 MM DRILL BIT WITH QUICK COUPLING X 5" FOR 4.5 MM CORTICAL & 6.5 MM CANCELLOUS SCREW	2
➤ BO1002INS5002	Ø 3.5 / 4.3 MM STEP DRILL BIT WITH QUICK COUPLING X 9" X 5.0 / 6.5 LOCKING CANCELLOUS SCREW	2
➤ BO1002INS5003	Ø 4.3 MM DRILL BIT WITH QUICK COUPLING X 9" FOR 5.0 MM LOCKING SCREW	2
➤ BO1002INS5004	Ø 4.5 MM DRILL BIT WITH QUICK COUPLING X 5" FOR 6.5 MM CANCELLOUS SCREW GLIDING HOLE	2
➤ BO1002INS5005	Ø 3.5 / 6.5 MM LOCKING DRILL SLEEVE FOR 6.5 MM LOCKING CANCELLOUS SCREW	2
➤ BO1002INS5006	Ø 4.3 / 5.0 MM LOCKING DRILL SLEEVE FOR 5.0 MM LOCKING & 5.0 MM LOCKING CANCELLOUS SCREW	2
➤ BO1002INS5007	Ø 3.2 / 4.5 MM DRILL & TAP SLEEVE COMBINED FOR 4.5 MM CORTICAL SCREW	1
➤ BO1002INS5008	Ø 3.2 / 6.5 MM DRILL & TAP SLEEVE COMBINED FOR 6.5 MM CANCELLOUS SCREW	1
➤ BO1002INS5009	Ø 4.5 MM BONE TAP WITH QUICK COUPLING FOR 4.5 MM CORTICAL SCREW	1
➤ BO1002INS5010	Ø 6.5 MM BONE TAP WITH QUICK COUPLING FOR 6.5 MM CANCELLOUS SCREW	1



4.5 / 5.0 MM LOCKING - INSTRUMENT SET

REF. NO	COMMODITY NAME	QTY
➤ BO1002INS5011	T-HANDLE WITH QUICK COUPLING	1
➤ BO1002INS5012	Ø 4.5 MM SCREW DRIVER WITH FIBRE HANDLE X 9"	1
➤ BO1002INS5013	Ø 1.5 MM THREADED GUIDE WIRE X 150 MM	2
➤ BO1002INS5014	Ø 2.0 MM THREADED GUIDE WIRE X 225 MM	2
➤ BO1002INS5015	Ø 4.5 / 6.5 MM DEPTH GAUGE	1
➤ BO1002INS5016	4.5 MM TORQUE LIMITING SCREW DRIVER X 9"	1
➤ BO1002INS5017	Ø 4.5 / 5.0 MM LOCKING INSTRUMENTS BOX WITH IMPLANT & SCREW TRAY	1



BIOMEDORTHO
define excellence



WORKS UNIT



Plot No 1 & 24
Raghuvver Industrial Estate,
Opp Jay Chemicals,
Odhav Road,
Ahmedabad,
Gujarat - 382415.



CORPORATE OFFICE



A/242, 2nd Floor,
Siddharth Annexe - 3,
Vemali - 391740,
Vadodara,
Gujarat - India.



FOR ENQUIRY

✉ info@biomedortho.com
✉ exports@biomedortho.com
☎ (+91) 74330 12333



www.biomedortho.co.in